

**PAJUNK®**

**NerveGuard**  
Automatic system  
for injection pressure limitation



Peripheral nerve blocks

MADE IN GERMANY


## Causes and approaches

# Avoiding nerve damage during peripheral nerve blocks

Ultrasound-guided localisation of peripheral nerves and the associated real-time visualisation provide crucial benefits in regional anaesthesia.<sup>1</sup> Nonetheless, it is clear that this does not reduce the incidence of permanent nerve damage.<sup>2</sup> Even in combination with nerve stimulation, intrafascicular injections cannot be ruled out.<sup>1,3,4,5</sup> The causes include incorrect positioning of the cannula along with exceeding critical injection pressures, as the following explanations show.


### Causes of the development of nerve damage

#### 1. Intrafascicular injections



The higher the injection pressure is during injections in the intrafascicular space, the more severe and longer lasting the resultant paresis is.<sup>6</sup>

If the pressure during intrafascicular injections exceeds a critical threshold of 15 psi, such injections may demonstrably lead to severe long-term neurological complications.<sup>6,7,8</sup> (Kapur<sup>6</sup>: > 20 psi, Hadzic<sup>7</sup>: > 25 psi, Hasanbegovic<sup>8</sup>: > 15.9 psi)



Intrafascicular injection

**Effects:** In the case of several hours lasting intrafascicular injections at high pressure, the microvascular blood supply of the nerve is severely restricted, which can lead to degeneration of nerve structures.<sup>7,12</sup>

→ *Limiting the injection pressure avoids nerve damage.*

1 Choquet, Capdevila, Ultrasound-guided nerve blocks ..., 2012 May; 114(5): 929–930

2 Neil et al., The Second American Society of Regional Anesthesia ..., 2016 March–April; 41(2): 183

3 Robards et al., Intraneural injection with low-current stimulation ..., 2009 Aug; 109(2): 673–677

4 Vassiliou et al., Risk evaluation for needle-nerve contact related ..., 2016 Mar; 60(3): 400–406


5 Sites et al., Characterizing novice behavior associated with ..., 2007 Mar–Apr; 32(2): 107–115

6 Kapur et al., Neurologic and histologic outcome ..., 2007 Jan; 51(1): 101–107

7 Hadzic et al., Combination of intraneural injection ..., 2004 September–October; 29(5): 417–423

8 Hasanbegovic et al., Effects of intraneural and perineural injection ..., 2013; 3(3): 248

## 2. Cannula-nerve contact




Incidence of paresthesias, motor evoked responses as well as exceeding the opening pressure threshold in the case of three different cannula positions.<sup>10</sup>

**Effects:** Direct cannula-nerve contact can also lead to damage of the neural structures with subsequent transient or permanent neurological impairment.<sup>10,11</sup> Localisation control using ultrasound and/or nerve stimulation may not in all cases reliably indicate direct cannula-nerve contact.

Paresthesia in and of itself is not regularly observed.<sup>10</sup> Avoiding direct cannula-nerve contact minimises the risk of damage to the nerve wall.<sup>11</sup>

→ *A reliable indicator of direct cannula-nerve contact is a high opening pressure.*




Direct cannula-nerve contact at high opening pressure.

## 3. Cannula-fascia contact

**Effects:** Misinjections into the wrong tissue layers can also be the cause of anaesthesia failure. For example, a high opening pressure can indicate that the tip of the cannula is closed by an in front fascia.<sup>9</sup>

→ *One approach to avoid misinjections into wrong nerve tissue structures is to limit the opening pressure.*



Blocking of the cannula opening by the fascia lata at high opening pressure.

9 Gadsden et al., High Opening Injection Pressure is Associated With Needle-Nerve ..., 2016 Jan–Feb; 41(1): 50–55

10 Gadsden et al., Opening injection pressure consistently detects needle-nerve contact ..., 2014 May; 120(5): 1246–1253

11 Steinfeldt et al., Histological consequences of needle-nerve contact following nerve stimulation ..., 2011; Article ID 591851: 0–9

12 Lundborg et al., Nerve compression injury and increased endoneurial fluid pressure ..., 1983 Dec; 46(12): 1119–1124

## NerveGuard


# The automatic injection pressure limiter

With the automatic injection pressure limiter NerveGuard, PAJUNK® addresses the issue of preventive position control. If the opening or injection pressure reaches the specified threshold, the administration of anaesthetic is automatically stopped. In this way, incorrect positioning of the cannula can be detected and corrected immediately and in addition an injection performed using too much pressure can be avoided as far as possible.

- If the pressure rises too high, the injection blocks automatically
- No monitoring and no visual control necessary
- Suitable for single shot and continuous peripheral nerve blocks



NerveGuard  
Injection pressure limiter



If the system pressure reaches the preset limit value, the NerveGuard valve closes automatically and the anaesthetic supply is blocked.




PAJUNK® cannula


The user can concentrate on the puncture and the ultrasound device.  
→ No additional "eye contact" with the NerveGuard is necessary.

## Preventive position control through NerveGuard

NerveGuard offers additional support in the localisation of the cannula tip and prevents nerve damage, together with the automatic pressure limitation.




**Cannula-fascia contact  
during interscalene blocks**  
Opening pressure  $\geq$  Limit value  
Indicates occlusion of the tip of the  
cannula as a result of fascia in front of it.<sup>9</sup>



**Cannula-nerve contact  
during femoral nerve blocks**  
Opening pressure  $\geq$  Limit value  
Indication of direct cannula-nerve  
contact.<sup>10</sup>


### Pressure measurement reaches limit value

Valve of the NerveGuard closes automatically.




- Administration of the anaesthetic is noticeably blocked.
- As soon as the pressure drops when the cannula is withdrawn, the valve reopens.

### Position correction through withdrawal/repositioning of the cannula



→ The tip of the cannula pierces  
the fascia and is now in the optimal  
position for an injection.



→ Avoiding an injection when there is direct  
cannula-nerve contact minimises the risk of  
damage to the nerve wall.<sup>11</sup>

→ *The injection is successful – the anaesthetic spreads around the nerve.*


9 Gadsden et al., High Opening Injection Pressure is Associated With Needle-Nerve ..., 2016 Jan–Feb; 41(1): 50–55

10 Gadsden et al., Opening injection pressure consistently detects needle-nerve contact ..., 2014 May; 120(5): 1246–1253

11 Steinfeldt et al., Histological consequences of needle-nerve contact following nerve stimulation ..., 2011; Article ID 591851: 0–9

*All-round solution*


## ***The connection between preventive position control and precise localisation***



- Detects incorrect cannula positioning (fascia or nerve contact)
- Prevents intrafascicular injections in cases of excessive pressure
- Automatically blocks injections when injection pressure exceeds the limit value
- Provides important information regarding position correction
- Avoids nerve damage and increases patient safety

Automatic injection pressure limitation with NerveGuard

*"Things just got a whole lot smarter."*  
**NerveGuard – just to be safe!**



High-precision stimulation using NanoLine thin layer technology



Excellent echogenicity with ultrasound cannulas  
(Cornerstone reflectors)



Preventive position control to avoid nerve damage with NerveGuard

*Three innovative techniques to protect the nerves during injections:  
NanoLine thin layer technology | Cornerstone reflectors | NerveGuard =*

 NERVE  
PROTECT<sup>3</sup>

# *Studies*

- Choquet O., Capdevila X. Ultrasound-guided nerve blocks: the real position of the needle should be defined, *Anaesth. Analg.* 2012 May; 114(5): 929–930
- Gadsden J., Choi J.J., Lin E., Robinson A. Opening injection pressure consistently detects needle-nerve contact during ultrasound-guided interscalene brachial plexus block, *Anesthesiology* 2014 May; 120(5): 1246–1253
- Gadsden J., Latmore M., Levine D.M., Robinson A. High Opening Injection Pressure is Associated With Needle-Nerve and Needle-Fascia Contact During Femoral Nerve Block, *Reg. Anesth. Pain Med.* 2016 Jan–Feb; 41(1): 50–55
- Helen L., O'Donnell B. D., Moore E. Nerve localization techniques for peripheral nerve block and possible future directions, *Acta Anaesth. Scandi.* 59 2015; 962–974
- Kapur E., Vuckovic I., Dilberovic F., Zacicagic A., Cosovic E., Divanovic K.A., Mornjakovic Z., Babic M., Borgeat A., Thys D.M., Hadzic A. Neurologic and histologic outcome after intraneural injections of lidocaine in canine sciatic nerves, *Acta. Anaesthesiol. Scand.* 2007 Jan; 51(1): 101–107
- Lundborg G., Myers R., Powell, H. Nerve compression injury and increased endoneurial fluid pressure: a „miniature compartment syndrome”, *J. Neurol. Neurosurg. Psychiatry* 1983 Dec; 46(12): 1119–1124
- Neil J.M., Brull R., Horn J.L., Liu S.S., McCartney C.J., Perlas A., Salinas F.V., Tsui B.C. The Second American Society of Regional Anesthesia and Pain Medicine Evidence-Based Medicine Assessment of Ultrasound-Guided Regional Anesthesia, *Reg. Anaesth. Pain Med.* 2016 March–April; 41(2): 181–194
- Patil J., Ankireddy H., Wilkes A., Williams D., Lim M. An improvised pressure gauge for regional nerve blockade/anesthesia injections: an initial study, *J. Clin. Monit. Comput.* 2015 Dec; 29(6): 673–679
- Patil J.J., Ford S., Egeler C., Williams D.J. The effect of needle dimensions and infusion rates on injection pressures in regional anaesthesia needles: a bench-top study, *Anaesth.* 2015; 70: 183–189
- Robards C., Hadzic A., Somasundaram L., Iwata T., Gadsden J., Xu D., Sala-Blanch X. Intraneural injection with low-current stimulation during popliteal sciatic nerve block, *Anesth.* 2009 Aug; 109(2): 673–677
- Ross S., Edwards K., McFadden K., Bigeleisen P.E., Orebaugh S.L. Pressures of Injection in a Cadaver Model of Peripheral Nerve Blockade, *J. Anesth. Clin. Res.* 2014; 5: 10
- Sites B.D., Spence B.C., Gallagher J.D., Wiley C.W., Bertrand M.L., Blike G.T. Characterizing novice behavior associated with learning ultrasound-guided peripheral regional anesthesia, *Reg. Anesth. Pain Med.* 2007 Mar–Apr; 32(2): 107–115
- Steinfeldt T., Graf J., Schneider J., Nimphius W., Weihe E., Borgeat A., Wulf H., Wiesmann T. Histological consequences of needle-nerve contact following nerve stimulation in a pig model, *Anesth. Research and Practice* 2011; Article ID 591851: 0–9
- Vassiliou T., Müller H.H., Limberg S., De Andres J., Steinfeldt T., Wiesmann T. Risk evaluation for needle-nerve contact related to electrical nerve stimulation in a porcine model, *Acta Anaesthesiol. Scand.* 2016 Mar; 60(3): 400–406

# NerveGuard

## All information at a glance

### NerveGuard



- NerveGuard injection pressure limiter, single, sterile. For single shot applications in diameters of **20 G/21 G/22 G**

Art. no.	PU
001151-38M	10




- NerveGuard injection pressure limiter, single, sterile. For single shot applications in diameters of **24 G/25 G**

Art. no.	PU
001151-38N	10

### Compatible with PAJUNK® cannulas

Product	Size	Art. no.	PU	Product	Size	Art. no.	PU
<b>UniPlex NanoLine</b>							
Facet grinding	22 G x 40 mm	001156-70	10	SonoPlex	24 G x 40 mm	001185-30G	10
	22 G x 50 mm	001156-74	10		22 G x 50mm	001185-31G	10
	22 G x 80 mm	001156-71	10		22 G x 70mm	001185-31H	10
	22 G x 100 mm	001156-84	10		22 G x 90mm	001185-31J	10
	22 G x 120 mm	001156-82	10	SonoBlock	25 G x 50 mm	001180-81	10
	21 G x 80 mm	001156-88	10	Facet grinding	24 G x 40 mm	001180-78	10
	21 G x 100 mm	001156-77	10		24 G x 50 mm	001180-85	10
	20 G x 120 mm	001156-72	10		22 G x 40 mm	001180-70	10
	20 G x 150 mm	001156-76	10		22 G x 50 mm	001180-74	10
<b>UniPlex NanoLine</b>							
SPROTTE®	22 G x 40 mm	001156-31G	10		22 G x 80 mm	001180-71	10
	22 G x 70 mm	001156-31H	10		21 G x 100 mm	001180-77	10
	22 G x 90 mm	001156-31J	10		20 G x 120mm	001180-72	10
	22 G x 150 mm	001156-28L	10	SonoPlex	24 G x 40mm	001180-30G	10
<b>SonoPlex</b>							
Facet grinding	24 G x 25 mm	001185-75	10		22 G x 50 mm	001180-31G	10
	24 G x 40 mm	001185-78	10		22 G x 70mm	001180-31H	10
	22 G x 40 mm	001185-70	10		22 G x 90mm	001180-31J	10
	22 G x 50 mm	001185-74	10	SonoTAP	24 G x 40mm	1185-3Y040	10
	22 G x 80 mm	001185-71	10	Facet grinding	24 G x 50 mm	1185-3Y050	10
	21 G x 80 mm	001185-88	10		22 G x 50 mm	1185-3E050	10
	21 G x 100 mm	001185-77	10		22 G x 80 mm	1185-3E080	10
	20 G x 120 mm	001185-72	10		21 G x 110 mm	1185-3F110	10
	20 G x 150 mm	001185-76	10		21 G x 150 mm	1185-3F150	10



**PAJUNK® GmbH**  
Medizintechnologie  
Karl-Hall-Strasse 1  
D-78187 Geisingen / Germany  
Phone +49 (0) 77 04/92 91-0  
Telefax +49 (0) 77 04/92 91-6 00  
[www.pajunk.com](http://www.pajunk.com)

**PAJUNK® Medical Produkte GmbH**  
D.A.CH • BeNeLux  
Karl-Hall-Strasse 1  
D-78187 Geisingen / Germany  
Phone +49 (0) 77 04/80 08-0  
Telefax +49 (0) 77 04/80 08-150  
[www.pajunk.com](http://www.pajunk.com)